

MSCmedium

We can provide a culture medium suitable for your specific purpose.
Please feel free to contact us for a consultation.

Switch to
Animal origin
free



Enhance
cytokine
production

Inhibit cell
aggregation

Xeno-Free Medium

2023~
No human-
derived
components

Animal- Origin-Free Medium

2023~
No animal-
derived
components

Mammalian MSC Medium

2024~
For the proliferation
of canine-derived
MSCs

Cytokine Production Medium

2024~
Successfully enhanced the cells'
cytokine production

Blood Coagulation Inhibitory Medium

2024~
Inhibits the expression of CD142,
a cellular procoagulant factor

Medium for 3D Culture

2024~
Successfully suppressed cell
aggregation

MSC Establishment Medium

2025~
Optimized for MSC establishment


Since developing
serum-free medium
for MSCs in 2022, we
have continued
development tailored
to user needs.



Case 1: Inhibition of Cell Enlargement by Xeno-Free Medium


While hMSCs are observed to have a tendency toward cell hypertrophy upon repeated passaging, hMSCs cultured in our Xeno-Free Medium showed suppressed hypertrophy and maintained a more undifferentiated state compared to cells cultured in competitor media or FBS + GF medium.

*FSC (Forward scatter) = Indicates cell size.



Case 2: Animal-Origin-Free Medium


Out of more than 20 candidate media, two Animal-Origin-Free (AOF) media were able to support proliferation for 4 passages or more, and the surface markers showed high values compared to the control medium.



Case 3: Culture of Canine MSC Lines


We enabled the culture of canine amniotic membrane MSC lines in serum-free medium, previously requiring FBS serum medium.

Competitor serum free Medium




A high number of dead cells are observed in suspension.

our serum free Medium



Showed good proliferation.

Case 4: Evaluation of Cytokine Production



Umbilical cord-derived MSCs were cultured in T25 flasks using a competitor medium, our medium ①, and our medium ② for four days. The culture supernatant was then collected, and the cytokine concentration was measured using ELISA.

Compared to the competitor medium, our medium ① and our medium ② showed higher cytokine concentrations in the culture supernatant, and also tended to show higher production per cell count. Regarding the production of Cytokine A, our medium ② showed a higher tendency than our medium ①, while for Cytokine B, our medium ① showed a higher tendency.

Case 5 : Shortening the Manufacturing Timeline via Refined MSC Establishment and Culture Protocol from Adipose Tissue

The conventional method for recovering over 100 million MSCs, starting from seeding a small amount of adipose tissue, typically takes 4 to 7 weeks. However, by using our improved establishment and expansion protocol and medium, it is now possible to culture them within 2 to 3 weeks (only 2 passages) without the use of special equipment.




adipose tissue



Establishment Process	Expansion Process
Establishment Medium 50 ~ 100 mL	Proliferation Medium 500 ~ 700 mL
7 ~ 14 days	5 ~ 10 days

Our MSC medium supported proliferation up to 1×10^9 cells after 16 days of culture.




Microscopic Comparison of Cultured MSCs Established from 0.2 g of Adipose Tissue: Competitor Medium vs. Our MSC Medium



Comparison of Cell Counts of MSCs Established from 0.2 g of Adipose Tissue in Competitor Medium and Our MSC Medium


Case 6: Development of a Culture Supplement to Suppress Tissue Factor (TF) Expression Related to the Prothrombotic Potential of MSCs

Tissue factor (TF; CD142), one of the extrinsic blood coagulation factors, is considered a contributing factor to the prothrombotic potential of MSCs. We screened over 3,000 small molecule compounds from our library to identify those effective in inhibiting TF. This resulted in the development of a culture supplement capable of suppressing TF expression even in FBS-containing medium. We performed the following three evaluations on Adipose-Derived MSCs (ADSCs) cultured in DMEM/10%FBS supplemented with 6 hit supplements (No. 6, 13, 15, 19, 24, 37) identified during screening:



Supplements No. 19, 24, and 37 significantly suppressed TF expression.


FCM Evaluation of Cell Surface TF Quantity



The amount of TF on the ADSC surface was decreased in the supplement-added group.

Clotting Time Measurement of ADSCs and Human Plasma (Clotting Assay)

ADSC cell suspension (1.5×10^5 cells/mL, containing citrate) was mixed with an equal volume of human blood plasma. The mixture was set in an automatic coagulation measuring device (CA-104, Sysmex), and the clotting time (seconds) after the addition of calcium chloride (recalcification) was measured.



CA-104, Sysmex

ADSCs cultured with Supplements No. 19 and 37 showed a significant extension of the coagulation time when mixed with human plasma, indicating a reduction in procoagulant potential.